



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6700



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.










MD™ FOAMER w/ACTUATOR (110 mL, 0.7 mL Metered Dose)

Plastic, White w/Clear View Window

DESCRIPTION: The portfolio of MD™ Dispensers are engineered for accuracy with a medical structure, yet convey a cosmetic appeal. The MD™ Dispensers will maintain the integrity of a formulation, while allowing patients to dose medications easily and accurately.

The MD™ foamer comes with a clear view window to inspect the content and its level. It also allows the patient to dispense formulations in a convenient and economical manner, while achieving the same result as the regular liquid pumps.

FEATURES:

-  Airless dispensing
-  Metered dosing
-  Excellent evacuation
-  No priming required
-  UV-resistant
-  Hygienic use – no bacterial contamination
-  Latex-Free
-  BPA-Free
-  PVC-Free

TECHNICAL SPECIFICATIONS:

Overall Dimension, Ø x L (mm / in)	47.0 x 156.5 / 1.85 x 6.16
Volume (mL / oz)	110.0 / 3.7
Dosage / stroke (mL / oz)	0.70 / 0.02