



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING






PRODUCT NO.: 6712



Test Kit, Personal Aseptic Technique Test ,PATT (Sterile) **5 units/case**

DESCRIPTION: This kit is easy to use for testing up to 5 individuals.

FEATURES:

-  20mL vials of TSB media
-  100mL partially filled mini-bags for added complexity
-  Log sheet included
-  Gummed labels included
-  Product conforms to USP standards

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM