



PRODUCT NO.: 6738



EZ SERIES AUTOCLAVE, FULLY AUTOMATIC (28.5 L), TUTTNAUER

Model EZ11Plus, Ø 11" × L 19.8" Chamber, 5 Trays, Cassette Capacity: 4 Full & 4 Half

DESCRIPTION: The EZ11Plus fully automatic autoclave helps practitioners meet today's challenging workloads and regulatory requirements by using steam as the sterilizing agent to completely sterilize surgical and non-surgical instruments. Steam sterilization of medical instruments is safe and effective. It is also designed to meet the most current sterilization standards such as ANSI / AAMI ST55

Comes with: Handle for the stainless steel trays, Cassette tray rack, Pouch rack

FEATURES:

- ✔ The newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing
- ✔ Innovative EZView Multi-Color Display offers real time cycle information that is easy to read, even from a distance. The EZView displays cycle parameters, temperature, and pressure in the chamber. During a cycle each stage is displayed with an active color wheel and text
- ✔ Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high efficiency air pump circulates filtered air through the chamber for fast drying. 0.2 µ HEPA filter (replaceable)
- ✔ 11" diameter × 19.8" deep chamber and flexible loading options. The exceptionally deep 19.8-inch chamber works in a variety of loading options to suit any office:
 - ✔ Use the 5 full sized trays for packs/pouches
 - ✔ Use the tray rack (included) to accommodate 4 full plus 4 half sized cassettes
 - ✔ Use the pouch rack(included) to encourage proper loading
- ✔ Front or top water fill options: Accessible EZFill front or top water fill option to enable easier filling of the reservoir
- ✔ USB Port enables users to easily transfer cycle data via a flash drive to a PC
- ✔ Network Port allows for connection to a local network to provide remote monitoring of the system when used with optional RPCR Software. Network port is located in the rear of the unit
- ✔ Dynamic Air Removal Technology to effectively remove air from the chamber
- ✔ Designed to meet the most current sterilization standards ANSI/AAMI ST55
- ✔ Optional printer documents a detailed history of each cycle
- ✔ Multiple languages

TECHNICAL SPECIFICATION:

Chamber Dimensions, Ø x L (cm / in)	28.0 x 50.4 / 11.0 x 19.8
Overall Dimensions, W x D x H (cm / in)	53.0 x 63.0 x 44.0 / 20.9 x 24.8 x 17.3
Volume (L / gal)	28.5 / 7.5
Standard Cassette Capacity	4 Full & 4 Half (loaded vertically)
Tray Dimensions, W x D x H (cm / in)	16.9 x 41.3 x 1.6 / 6.7 x 16.3 x 0.6
No. of Trays	5
Voltage / Frequency	120 V / 50 / 60 Hz
Power / Current	1400 W / 12 A
Shipping Weight (kg / lb)	56 / 122

ADDITIONAL INFORMATION:

- ✔ A minimum tabletop depth of 22" is required
- ✔ The standard full size cassette is 8" x 11" and the standard half size cassette is 8" x 5.5"

WARRANTY:

- ✔ Parts and labor: Two (2) years
- ✔ Chamber: Ten (10) years