



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING



**PRODUCT NO.: 6779**



## **Mold, Suppository Aluminum 2g (100 Cavity)**

**DESCRIPTION:** Aluminum suppository molds produce 2 gm suppositories in a universal shape for rectal or vaginal use. Quick and easy assembly - no tools required. Threaded posts and wingnuts secure the mold in precise alignment.

### **FEATURES:**

-  Dimensions, L x W x H (in): 11.75 x 3 x 1.5
-  Weight (lbs): 4

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)