




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO. 6805**



**TENSIONING WRENCH (Automatic Inversina 2L)**

For Inversina Automatic Tumbler Mixer 2L

**DESCRIPTION:** Tensioning wrench is needed for tightening the vessel with different sizes in the Inversina cage securely before operation.

**ADDITIONAL INFORMATION:**

**APPLICABLE PRODUCTS:**

-  6810-01: Automatic Inversina Tumbler Mixer