



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO. 6806



ELECTRIC FUSES (Automatic Inversina 2L) For Inversina Automatic Tumbler Mixer 2L

DESCRIPTION: Electric fuses in the automatic Inversina provide overcurrent protection. They interrupt excessive current and prevent further damage to the entire system.

ADDITIONAL INFORMATION:

APPLICABLE PRODUCTS:

-  6810-01: Automatic Inversina Tumbler Mixer