



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

## PRODUCT NO. 6807



### **CARBON BRUSH (Automatic Inversina 2L)** For Inversina Automatic Tumbler Mixer 2L

**DESCRIPTION:** This sliding unit comes with Inversina Automatic Tumbler Mixer 2L and is used to transmit electrical current to the rotating part.

#### **ADDITIONAL INFORMATION:**

#### **APPLICABLE PRODUCTS:**

-  6810-01: Automatic Inversina Tumbler Mixer