




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO. 6808



TENSION RINGS (Automatic Inversina 2L)

For Inversina Automatic Tumbler Mixer 2L

DESCRIPTION: Tension rings in automatic Inversina are connected to rotating ring which need to be tightened to fit a variety of mixing containers for example to secure a wide-neck flasks in the cage.

ADDITIONAL INFORMATION:

APPLICABLE PRODUCTS:

-  6810-01: Automatic Inversina Tumbler Mixer