



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6851



SUPERNUOVA+ STIRRING HOT PLATE, THERMO SCIENTIFIC

7.25" x 7.25", Digital

DESCRIPTION: This digital stirring hotplate offers precise stirring control, safety and temperature performance for routine protocols.

FEATURES:

- ✔ Slow-speed stirring, consistent speed control, and strong magnetic coupling
- ✔ HOT TOP warning system to protect from accidental burns with prominent display when heating surface is above 50 °C (122 °F)
- ✔ User-adjustable over-temperature protection
- ✔ Easy-to-read LED display for heating and stirring
- ✔ Raised display design protects electronics from spills
- ✔ Temperature calibration based on application needs
- ✔ Programmable for repeated experiments
- ✔ Timer for untended operation
- ✔ Stir bar, temperature probe, rod and clamps included

WARRANTY: 3 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	20.8 x 33.0 x 9.7 / 8.2 x 13.0 x 3.8
Top Plate Dimensions, W x D (cm / in)	18.4 x 18.4 / 7.25 x 7.25
Weight (Kg / lb)	5.6 / 12.3
Top Plate Material	Ceramic
Temperature Range (°C / °F)	30 – 450 (86 – 842)
Stirring Range (RPM)	50 – 1500
Conformity	EN 61000-3-2 EN 61000-3-3 EN 61326-1 EN 61010-1 EN 61010-2-010 EN 61010-2-051 CE (230 V Models)