



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6852



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

CIMAREC+ STIRRING HOT PLATE, THERMO SCIENTIFIC

4.25" x 4.25", Digital

DESCRIPTION: This digital stirring hotplate offers precise stirring control, safety and temperature performance for routine protocols.

FEATURES:

- ✔ Slow-speed stirring, consistent speed control, and strong magnetic coupling
- ✔ HOT TOP warning system to protect from accidental burns with prominent display when heating surface is above 50 °C (122 °F)
- ✔ Automatic over-temperature protection
- ✔ Easy-to-read LED display for heating and stirring
- ✔ Raised display design protects electronics from spills
- ✔ Stir bar included

WARRANTY: 3 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	13.0 x 25.9 x 9.8 / 5.1 x 10.2 x 3.9
Top Plate Dimensions, W x D (cm / in)	10.8 x 10.8 / 4.25 x 4.25
Weight (Kg / lb)	2.2 / 5.0
Top Plate Material	Ceramic
Temperature Range (°C / °F)	30 – 540 (86 – 1004)
Stirring Range (RPM)	50 – 1500
Conformity	EN 61000-3-2 EN 61000-3-3 EN 61326-1 EN 61010-1 EN 61010-2-010 EN 61010-2-051