



**PRODUCT NO.: 6854**



## PRINTER (EPSON)

For SteriSure Dry Heat Sterilizers

**DESCRIPTION:** Epson's TM-U220 printers are compact, reliable and fast, and offer all the easy-to-use features important to a wide variety of industries.

### FEATURES:

- Driver Support: Windows, Linux, OPOS, SDK for iOS and Android
- Easy to operate
- Faster printer speed
- Drop-in paper loading
- Two-color black and red printing
- Impactful logo in black
- Flexible Paper widths

### TECHNICAL SPECIFICATIONS:

<b>Print method</b>	9-pin, serial impact dot matrix
<b>Print font</b>	Font: 7 x 9/9 x 9 Column capacity (columns): 40/42 or 33/35 Character size: 1.2(W) x 3.1(H)/1.6(W) x 3.1(H)mm Character set: 95 Alphanumeric, 37 International, 128 x 12 Graphic Traditional/Simplified Chinese, Thai, Japanese, Korean Characters per inch: 17.8/16 cpi or 14.5/13.3
<b>Interfaces</b>	RS-232, Bi-directional parallel, Connect-It™ family Dealer Option: USB, 10 Base-T I/F
<b>Data buffer</b>	4 KB or 40 bytes
<b>Print speed</b>	4.7 lps (at 40 columns, 16 cpi), 6.0 lps (at 30 columns, 16 cpi)
<b>Paper dimensions</b>	57.5±0.5mm, 69.5±0.5mm, 76±0.5mm (W), x dia. 83.0
<b>Paper thickness</b>	0.06 to 0.085
<b>Copy capability</b>	One original and one copy
<b>Inked ribbon</b>	ERC-38 (Purple, Black, Black/Red)
<b>Ribbon life</b>	Purple: 4 x 106 characters, Black: 3 x 106 characters B/R: 15 x 105 characters (Black), 75 x 104 characters (Red)
<b>Power</b>	AC adapter (included)

<b>Power consumption</b>	24 V
<b>D.K.D. function</b>	2 drivers
<b>Reliability</b>	MTBF 18 x 104 hours MCBF 18 x 106 lines
<b>Overall dimensions (mm)</b>	Type A: 160 (W) x 286 (D) x 157.5 (H) Type B: 160 (W) x 248 (D) x 138.5 (H) Type D: 160 (W) x 248 (D) x 138.5 (H)
<b>Mass (Approx.)</b>	Type A: 2.7 Kg Type B: 2.5 Kg Type D: 2.3 Kg
<b>EMC standards</b>	VCCI class A, FCC class A, CE marking, AS/NZS 3548 class B
<b>Power supply</b>	Exclusive AC adapter (included)
<b>Factory option</b>	Near-end sensor
<b>Option</b>	Wall-hanging bracket WH-10
<b>Safety standards</b>	UL, CSA, TUV (EN60950)

**ADDITIONAL INFORMATION:**

**Warranty:** 2 years

**Applicable Products:**

-  7873: Dry Heat Sterilizer, 2100
-  7420: Dry Heat Sterilizer, 3100