



**PRODUCT NO.: 6859**



## **SETTLING PLATE, TSA (Sterile)**

25mL, 100 x 15mm, Double-Wrapped

**DESCRIPTION:** Its revolutionary design has been engineered to provide maximum security after collecting samples in critical areas. With a simple twist of the lid, the plate securely locks, thereby preventing accidental openings during handling and transport.

### **FEATURES:**

- Prevents contamination caused by accidental opening and provides additional security during the sampling process
- No need for unnecessary and risky elastic or adhesive wrapping to ensure the lid stays closed
- Lid remains securely attached to the base assuring the integrity of the samples and the performance of the media

### **ADDITIONAL INFORMATION:**

#### **Applicable Products:**

- 6861: Microbial Air Sampler, Sampl' Air Pro