



PRODUCT NO.: 6908



SNAP-ON STOPPER (13mm, Gray)

Chlorobutyl, Siliconized

DESCRIPTION: This design provides maximum sealing for vials with a blow-back glass finish, which is a small ridge on the interior of the vial's neck. The snap-on stopper locks into the blow-back ridge, providing a tight seal and keeping the stopper from backing out. This style is not recommended for vials with a straight-walled finish. The blow-back feature is typically found in vials from North American manufacturers.

FEATURES:

- Manufactured from the highest quality raw materials
- Components manufactured to yield exceptional dimensional stability
- Snap-on style provides maximum sealing for vials with a blow-back glass finish
- Free of latex, nitrosamines and 2-MCBT (2- mercaptobenzothiazole)
- Very well suited for ready-for-sterilization and ready- for-use applications
- Autoclavable

TECHNICAL SPECIFICATIONS:

- Durometer: 50
- Style: Snap-On
- Rubber Material: Gray Chlorobutyl, Siliconized