



PRODUCT NO.: 6911



SPILL TRAY, AIR SCIENCES (3ft)

Polypropylene, For Purair Basic Ductless FumeHood, Black

DESCRIPTION: This chemically resistant spill tray is used with the Purair Basic Ductless Fume Hood (6912-01) to contain spills and simplify cleaning procedures.

TECHNICAL SPECIFICATIONS:

- Dimensions, W x D x H (in / cm): 35 x 27 x 1 / 89 x 68.6 x 2.5
- Weight (lbs / kg): 15 / 6.8

ADDITIONAL INFORMATION:

Applicable Products:

- 6912-01: 3ft Purair Basic Ductless Fume Hood