



PRODUCT NO.: 6943







EVAPORATION FLASK, IKA

500mL, NS 29/32, Borosilicate Glass, Clear

DESCRIPTION: This evaporation flask can be used in conjunction with the IKA Grinder Drive Station to create a closed system when pulverizing your tablets into powder. The 500mL flask is made of high quality borosilicate glass.

ADDITIONAL INFORMATION:

Applicable Products:

-  6991: Grinder, Drive Station – REQUIRED
-  6990: Particle Sieve, 0.25mm – REQUIRED
-  6989: Impact Grinding Head – REQUIRED
-  6953: Clamp – REQUIRED