



PRODUCT NO.: 6976



PROMOLD-TLM, TORPAC

100 Cavities x 0.63mL, Hard Anodized Aluminum, For Making Troches & Lozenges, Black

DESCRIPTION: Torpac's Hard Anodized Aluminum ProMold-TLM is an ergonomic and effective solution for making troches and lozenges.

FEATURES:

- For making molded Lozenges & Troches
- 100 holes/cavities
- Volume: 0.63ml or 725mg Polyethylene Glycol 1450

TECHNICAL SPECIFICATIONS:

Overall Dimensions, W x D x H (in / cm)	12.5 x 6.0 x 0.6 / 31.6 x 15.2 x 1.45
Individual Mold Hold, Ø x H (in / mm)	0.5 x 0.2 / 12.7 x 5
Individual Troches/Lozenges, Ø x H (in / mm)	0.48 x 0.20 / 12.3 x 5.0
Materials of Construction	Hard Anodized Aluminum Inox 316 (Stainless Steel) Helicoils for screws Inox 316 (Stainless Steel) Screws

ADDITIONAL INFORMATION:

Cleaning Instructions:

- Clean with mild dishwashing liquid and hot water
- Dishwasher safe when cleaned with water only - no automatic dishwasher chemicals.