



PRODUCT NO.: 6990



GRINDER, PARTICLE SIEVE (IKA)

0.25mm, Stainless Steel, For IKA Grinder Drive Station

DESCRIPTION: Interchangeable sieve for insertion into the grinding heads ensures maximum particle size filtering. Use this product with drive station (6991-01) and grinding head (6989-01) to pulverize your tablets into powder.

FEATURES:

- Particle Sieve (0.25mm) allows to control particle size of final product

ADDITIONAL INFORMATION:

Applicable Products:

- 6991: Grinder, Drive Station – REQUIRED
- 6989: Impact Grinding Head – REQUIRED
- 6943: Evaporation Flask, 500mL – REQUIRED
- 6953: Clamp – REQUIRED