



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7056









BEAKER HOLDER

13"L x 5.5"Ø, Adjustable, Stainless Steel

DESCRIPTION: These sturdy, easily adjusted, multi-purpose clamps are ideal for holding large cylindrical glassware and similar objects. A large, flat thumbscrew/worm-drive permits maximum tightening or removing in minimum time. Clamp is constructed of corrosion-resistant stainless steel. Use this product in conjunction with 7584-01 and 8448-01 in order to hold a flask or beaker on a hot plate stirring unit.

FEATURES:

-  Material: Corrosion Resistant Stainless Steel
-  Grip Size (Min-Max): 2.5-5.5 in / 64-140 mm
-  Maximum Grip Size: 5.5" / 140 mm
-  Arm Length: 8.3" / 211 mm
-  Arm Diameter: 0.43" / 11 mm
-  Overall Length: 13.3" / 338 mm

ADDITIONAL INFORMATION:

Related Products:

-  7584-01: Hot Plate Support Rod
-  8448-01: Hook Connector