



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7065



4-PORT MANIFOLD, LABCONCO For FreeZone Freeze Dry System

DESCRIPTION: This sits on the freeze dryer and allows for connection, via the valves, of flasks using 1/2" and 3/4" adapters during the lyophilization process.

TECHNICAL SPECIFICATIONS:

| | | |
|---------------------------------|------------|---------------------------------------|
| Dimensions, W x D x H (cm / in) | | 21.8 x 21.8 x 28.4 / 8.6 x 8.6 x 11.2 |
| Weight (kg / lb) | | 2.5 / 5.5 |
| Material of Construction | Manifold | 304 Stainless Steel |
| | Valve | Neoprene |
| | Valve Knob | Plastic |