



PRODUCT NO.: 7080



MD PEN

6.5mL, 0.10mL Metered Dose, White w/Clear View Window

DESCRIPTION: The new portfolio of MD Dispensers was engineered for accuracy with a medical structure, yet designed to convey a cosmetic appeal. The MD Dispensers will maintain the integrity of a formulation, while allowing patients to easily and accurately dose medications. The first of its kind on the market, these dispensers feature UV Resistant view windows for added convenience.

FEATURES:

- Airless dispensing
- Slim, sleek, pen-type look
- Silicone tip for smooth application
- Excellent evacuation
- Minimal priming required
- UV-resistant
- Protecting from oxidization – extended product lifetime
- Hygienic use – no bacterial contamination
- Latex Free
- BPA Free
- PVC Free

TECHNICAL SPECIFICATIONS:

Fill Capacity (mL)	6.5
Dosage (mL)	0.10
Precision (% RSD)*	< 3.0
Materials of Construction	Piston: Polyethylene Inner Bottle: Polypropylene Outer Bottle: Polypropylene Cap: Polypropylene Nozzle: Polypropylene Tip: Silicone
Dimensions: L x Ø (in / cm)	5.5 X 0.8 / 13.9 x 2.0
Filling Method	3 Piece Assembly: 1) Top-Filled Pen 2) Snap-Cap Spring Pump 3) Outer Casing
Color	White w/Clear View Window
UV-Resistance	Bottle: Yes View Window: Yes


*Values established using MEDISCA's Cream Bases.

ADDITIONAL INFORMATION:

Filling Instructions:

MD Pen

TO FILL :



TO USE :

