



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7081**



## **STRAINER**

**2.5"L x 2.8"Ø, 40 Mesh, Stainless Steel Mesh w/Plastic Handle, White**

**DESCRIPTION:** Small strainer used for straining powders.

### **FEATURES:**

- ▣ Stainless Steel Mesh – 40 Mesh size
- ▣ Frame and handle made of durable plastic
- ▣ Features 2 "lips" on one side and handle on the other for conveniently resting on bowl, if needed
- ▣ Measures 2.5"/6cm diameter