



PRODUCT NO.: 7110-01



## BUBBLE POINT TEST APPARATUS (Nalgene)

For Nalgene Vacuum Filter Units, Cast Aluminum, Blue

**DESCRIPTION:** This apparatus permits convenient confirmation of membrane integrity and pore size of Thermo Scientific™ Nalgene™ MF75 Series Filter Unit or Bottle Top Filters directly after filter unit use without removing the membrane from the filter unit. Tests filter unit membranes without removing the membrane from the filter unit.

### FEATURES:

- Cast aluminum
- Tests the bubble point of 50, 75 and 90mm Nalgene filtration products
- Calibrate gauge and pressure source not included
- Accepts 0.25" (6mm) I.D. tubing
- Includes: Adapters and O-rings