



PRODUCT NO.: 7166-01



MIXING VESSEL (Inversina)

1.5L, PET, Clear w/White Screw Cap

DESCRIPTION: This 1.5L plastic mixing vessel can be used with the Inversina Manual Tumbler Mixer for mixing purposes. This product is not dishwasher safe. Dimensions (Ø x H): 4.7 x 6.5 in / 12 x 16.5 cm

ADDITIONAL INFORMATION:

Applicable Products:

-  7173-01: Inversina Manual Tumbler Mixer

Required Accessories:

-  7158-01: Band Clamps, 1.5L