



PRODUCT NO.: 7181-01



CHEW-A-TREAT® FISH MOLD

60 Cavities x 1.0mL, PVC

DESCRIPTION: MEDISCA's Chew-A-Treat® product line represents a simple, unique, and convenient avenue for dispensing chewable animal treats for your beloved pet. Now introducing the Chew-A-Treat® solution for your cat! This mold allows the manufacture of medicated biscuits with high palatability and appearance similar to veterinary cookies. Discover the solution veterinarians have been waiting for!

FEATURES:

- ✔ Individualized doses: 60 Cavities of 1.0mL
- ✔ Adapting the formulas for the discerning palate of cats
- ✔ Custom fish-shaped format for cats
- ✔ Includes a convenient tablet ejector
- ✔ Refrigeration no longer needed to removed treats from molds
- ✔ Instruction Manual Included

TECHNICAL SPECIFICATIONS:

Material	PVC
Cavity Plate Dimensions (mm)	235 x 5.5 x 258
Base Plate Dimensions (mm)	235 x 8 x 258
Extractor Dimensions (mm)	60 x 40 x 206
Side Clamps	25 x 45 x 257