



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7188-01







DRY BLOCK INCUBATOR

55–60°C, 15 x 11mm Cavities w/Vial Crusher, 110V, Aluminum

DESCRIPTION: This dry block incubator can easily be plugged in for monitoring the steam sterilization process with the self-contained biological indicators (7894-01).

FEATURES:

-  Temperature Range: 55-60°C / 131-140°F
-  Material: Aluminum
-  Cavities: 15 wells, 11mm each
-  Electrical Requirements: 110V