



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7201-01



AIR-TIGHT DAMPER (Labconco)





10"H x 10"Ø, Stainless Steel, For Labconco Purifier Biological Safety Cabinet

DESCRIPTION: Mounts atop the Canopy Connection to adjust exhaust airflow and closes for decontamination procedures. Dampers adjust the biosafety cabinet's exhaust airflow.

Applicable products:

-  8829-01
-  7588-01
-  7202-01
-  7200-01

FEATURES:

-  Manually adjustable
-  Type 304 Stainless Steel
-  10" (25.4cm) high with 10" (25.4cm) OD inlet and outlet
-  Weight: 13.0 lbs / 5.9 kg