



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7285**



## **pH METER POCKET PRO+, HACH (pH 0.0 to 14.0) w/Replaceable sensor**

**DESCRIPTION:** This pH tester is used to measure pH of general water samples.

Included in box: Tester (1), AAA Batteries (4), Lanyard (1), Single-Use 7.00 Calibration solution (1) and Quick start guide (1).

### **FEATURES:**

- ✓ Easy calibration steps and built-in diagnostics
- ✓ Easy to read LCD screen
- ✓ Comes with 4 easy-to-replace AAA batteries and a sample cup.

**WARRANTY:** 1 year for tester and 6 months for replacement sensor

### **TECHNICAL SPECIFICATIONS:**

Dimensions, W × D × H (mm / in)	37 × 30 × 170 / 1.45 × 1.18 0 × 6.69
Weight, with batteries (g / lb)	135 / 0.297
pH Range	0.0 – 14.0
pH Resolution	0.01
pH Accuracy	± 0.02
Calibration points	3 (Default); 2 (Custom)
Auto recognition calibration standards	USA: 4.01, 7.00, 10.01 pH; NIST: 4.01, 6.86, 9.18 pH
Power requirements	AAA alkaline batteries (4x)
Enclosure rating	IP67
Certifications	CE, FCC, Industry Canada, KC, RCM, China RoHS