



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7310-01**



## SIEVE

50 Mesh, 8" Diameter, 300 $\mu$ , Brass



**DESCRIPTION:** The sieve serves to control particle size distribution as well as the flow property of the blend. This is to ensure content uniformity of the solid dosage form.

### TECHNICAL SPECIFICATIONS:

Sieve Diameter (in / cm)	8 / 20.3
Sieve Material	Frame, Cloth: Brass
Sieve Height, H x D (in / cm)	2.625 x 2 / 6.7 x 5.1
U.S. Standard	No. 50
Metric Size ( $\mu$ , micron)	300

### TECHNICAL SPECIFICATIONS:

#### Applicable Products:

-  9652-01
-  9653-01