



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7315-01



HOT HAND PROTECTOR

Silicone Rubber, -57 to 260°C / -70 to 500°F, One Size Fits All, Red

DESCRIPTION: The Hot Hand® Protector mitt provides superior protection when handling hot or cold objects. The gripping surface is studded with suction cups for a secure no-slip grasp to safely remove the Medi-RDT Mold from the convection oven.

FEATURES:


- ✔ Use for handling hot and cold objects from -57°C (-70°F) to 260°C (500°F)
- ✔ Soft and compliant silicone rubber naturally conforms to hand and object gripping surfaces studded with small concave nubs for positive grip and additional protection against temperature extremes
- ✔ Easier and faster to use than glove or tongs for many applications
- ✔ Tab with hole for hanging hand protector near work station
- ✔ One size fits all hands, left or right 10 x 19cm (3.875 x 7.5 in.)

TECHNICAL SPECIFICATIONS:

Temperature Range (°C / °F)	-57 to 260 / -70 to 500
Material of Construction	Silicone Rubber w/small concave nubs for positive grip
Dimensions (in / cm)	3.875 x 7.5 / 10 x 19
Size	One size fits all, Left or Right

ADDITIONAL INFORMATION:

Related Products*:

-  2772-05: Medi-RDT Base, 100g
-  7219-01: Medi-RDT Blister Pack, 160 Doses
-  7219-02: Medi-RDT Blister Pack, 800 Doses
-  7314-01: Medi-RDT Mold, 96 Cavities x 750mg
-  7310-01, 9652-01: Sieve & Receiver Pan, 50 Mesh
-  9675-01: Powder Scraper

***Please note:** 7218-01 is the Medi-RDT Kit which will provide a complete package for RDT Molds

Required Accessories:

-  9011-01: Convection Oven, 70.8L / 2.5 cuft

Optional Accessories:

-  7311-01: Medi-RDT Multi-Dose Mold (180 Cavities: 60 x 75 mg, 60 x 150 mg, 60 x 200 mg)