



PRODUCT NO.: 7363



DOT MATRIX PRINTER W/CABLE, OHAUS

For AV, DV, EX, PA & SP* Pro Series Balances

DESCRIPTION: The all-new SF40A combines high performance and multiple functions to offer OHAUS quality and durability in a portable dot-matrix impact printer. The SF40A offers advanced functionality in an easy-to-use and set up printer at an exceptional value. This versatile printer is ideal for use in conjunction with OHAUS balances and scales in laboratory, industrial, and educational applications to ensure accurate documentation of weighing results.

Compatible with: Ohaus Discovery, Explorer, Pioneer, Scout Pro* & Adventurer Pro Series Balances.

*Required Interface Kit must be purchased alongside the appropriate balance. Please inquire if you hold a Scout Pro Balance.

FEATURES:

- Paper Roll, Ink Ribbon, Cable and Cable Adapter included
- Multiple Application Modes and Features Support the Data Reporting Process: SF40A is equipped with statistics and totalization modes that provide additional functionality for QA/QC and reporting processes. SF40A also offers a real-time clock with battery backup and date and time on printouts.
- Easy Setup And Configuration Eliminate Wasted Time: Don't let your printer inhibit the efficiency of your weighing process. OHAUS balances and scales improve productivity and increase throughput and SF40A speeds recordkeeping to keep your job on track and get it completed as proficiently as possible. PC configuration software makes connecting the printer between balance or scale and computer a fast and convenient process and also makes setting up date, time, and baud rate hassle-free. The SF40A combined with GLP/GMP data output can help meet traceability and compliance requirements.
- Multi-Language Printing Supports Universal Operation: SF40A supports printing in the language of the connected device in 13 different languages.

TECHNICAL SPECIFICATIONS:

Model	SF40A
Number of Characters per Line	40
Print Speed	Approximately 1.2 lines per second
Application Modes	Totalization, Statistics (Statistics printout includes number of samples, average value, standard deviation, minimum, maximum, min-max differential, and total sample weight sum)
Paper Roll (W x Dia.)	Standard paper 2.26 in. to 2.36 in. Ø max / 57.5 mm to 60 mm Ø max.
Power Requirements	AC adapter Input: 100 – 240 VAC, 50 - 60 Hz AC adapter output: 12 VDC, 2.5 A
Interface	RS232C
Ambient Operating Temperature	32 °F – 104 °F / 0 °C – 40 °C
Printer Dimensions (W x L x H)	5.5 x 8.1 x 2.9 in. / 140 x 205 x 75 mm
Net Weight	1.6 lb / 0.7 kg (Paper roll included)
Shipping Weight	1.8 lb. / 0.8 kg

ADDITIONAL INFORMATION:

Warranty: 1 year

Applicable Products:

- ☑ Ohaus Adventurer AX Series: 6942, 6874, 6941, 6940, 6939, 6938, 6937, 6936
- ☑ Ohaus Adventurer Pro Series (discontinued): 7280, 7212, 7783, 8656, 8938, 7481, 7341, 8204
- ☑ Ohaus Discovery Series: 8906
- ☑ Ohaus Explorer Series: 7760, 7810, 8016, 7968, 8014, 7276, 7793, 7853, 6951, 7555
- ☑ Ohaus NVL Series: 7215
- ☑ Ohaus Pioneer Series: 6878, 7768, 6906, 8194