



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7398



DISPERSING ELEMENT, IKA (10 – 1500 mL)

For T 25 Digital Dispersing Instrument

DESCRIPTION: This dispersing element is an accessory to the T 25 Dispersing instrument. It functions to properly homogenize and emulsify solutions and prepare good suspensions.

TECHNICAL SPECIFICATIONS:

Diameter, Stator / Rotor (mm)	18.0 / 12.7
Shaft Length (mm)	194
Immersion Depth (mm)	40 – 165
Working Range, Water (mL)	10 – 1500
Material of Construction	PTFE and AISI 316L
Maximum Allowable Speed (RPM)	25000
Maximum Temperature (°C)	180
Ultimate Fineness, Suspensions (µm)	10 – 50
Ultimate Fineness, Emulsions (µm)	1 – 10

WARRANTY: 2 years