



PRODUCT NO.: 7484



SERUM BOTTLE, NON-STERILE

30mL, 20mm, Glass Type 1, Pre-Cleaned & Depyrogenated, Amber

DESCRIPTION: The blowback feature in the neck and specially designed bottom radius, for added strength, are important for lyophilization applications. The clear vials are manufactured from low extractable borosilicate glass that conforms to ASTM Type I, Class A and USP Type I requirements. The amber vials are manufactured from low extractable borosilicate glass that conforms to USP Type I requirements. Vials are packaged in shrink-wrapped modules to greatly reduce particulate contamination.

FEATURES:

- Made from low extractable borosilicate glass
- Special design adds strength for freeze-drying applications
- Shrink-wrapped modules reduce particulate contamination

TECHNICAL SPECIFICATIONS:

Bottle Material	Glass
Bottle / Vial Style	Serum
Capacity (mL / oz)	30 / 1.01
Diameter (mm)	37
Height (mm)	62
Bottle Color	Amber
Mouth ID x OD (mm)	13 x 20
Glass Type	Type 1 Borosilicate
Finish Type	Crimp/Serum
Cap/Closure Size (mm)	20