



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7542



Balance Accessories Calibration Weight, 1.0g

DESCRIPTION: Can be used as a reference standard in calibrating other weights and is appropriate for calibrating high-precision analytical balances with a readability as low as 0.1 mg to 0.01 mg.

FEATURES:

- ✔ Tolerance: 0.034mg
- ✔ Style: Cylindrical with handling knob
- ✔ Material: 316 Stainless Steel
- ✔ Density: 7.95 g/cm³ at 20°C
- ✔ Dimensions (Weight Only): 0.184" x 0.417"
- ✔ Dimensions (Weight in Case): 2" x 2.5"
- ✔ Case: Supplied in high quality polycarbonate case with red insert

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM