



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING








PRODUCT NO.: 7543



Balance Accessories Calibration Weight, 0.01g

DESCRIPTION: Can be used as a reference standard in calibrating other weights and is appropriate for calibrating high-precision analytical balances with a readability as low as 0.1 mg to 0.01 mg.

FEATURES:

-  Tolerance: 0.010mg
-  Style: Flat with one end turned up for easy handling
-  Material: 304 Stainless Steel
-  Density: 7.95 g/cm³ at 20°C
-  Dimensions (Weight Only): 0.184" x 0.153" x 0.080"
-  Dimensions (Weight in Case): 2" x 2.5"
-  Case: Supplied in high quality polycarbonate case with red insert

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM