



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7708**









## MEGA PUMP AIRLESS DISPENSER

150mL, MACRO w/Cap, Compact, Bottom-Fill, 1.0mL Dosage, Clear

**DESCRIPTION:** The Mega Pumps are most suitable for creams, gels, pastes and lotions and ideal for products with natural and active ingredients. Dispensing is ideal for medium to high viscosity products. The Mega-pumps will provide optimal protection from exposure to air.

### FEATURES:

-  UV-resistant
-  Reliable dual-valve system
-  Dispenses in any position (360°)
-  Precise and repeatable dosage
-  100% plastic – No metal parts; Recyclable with mixed plastics
-  Each dispenser is 100% function tested

### TECHNICAL SPECIFICATIONS:

Dispensing Rate (mL)	1.0
Filling Method	Bottom-Fill
Color	Clear
Shape	Compact
Diameter (mm)	47
Height w/ cap (mm)	177

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)