



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING









PRODUCT NO.: 7919



APT-COM™ 3, BASIC-EDITION, BINDER

DESCRIPTION: The system offers features that are needed for tasks ranging from the simplest measurements to guideline-compliant work: monitoring of processes and documentation of process data. Documentation is automatically generated in electronic format and as a hard copy. This produces guideline-compliant documentation without extra effort. Easy to use for a broad range of applications.

FEATURES:

-  Remote control for connected simulation chamber
-  Graphic recording of measured values
-  Documentation of measured values in a tamper-proof file
-  Documentation when an adjustable tolerance limit is exceeded
-  Documentation of all system service with username and supporting documents
-  Blocking function (software lock) for unsupervised operation
-  Graphical/numerical program editor
-  Output of measured values in ASCII format on a printer or in a spreadsheet file for further processing, e.g. in spreadsheet programs, date range freely selectable