



**PRODUCT NO.: 8006**



## **PROFILLER 3700 CAPSULE FILLING SYSTEM, TORPAC (SIZE 1)**

### **300-Hole Capsule Machine**

**DESCRIPTION:** ProFiller 3700 is the new standard in 300-hole bench-top capsule filling systems. Experienced users like the higher output, fast size changes, easy to clean design & optional accessories. First time users love the ease-of-use & upgrade options.

Everything you need to fill one capsule size: Hand-Held Filler Base, Orienter Base and Locking Plate / Change Parts for one capsule size: Caps Tray, Body Sheet Set (A, B, C), Spacer Plate, Capsule Locker, Orienter Tray / 300-pin tamper / Powder Spreader / Powder Tray + Powder Tray Clamps / Illustrated Manual / Plastic storage box / Empty and Locked Capsule Bins / Set of empty capsule measuring bowls

#### **FEATURES:**

- Shipped fully assembled. Ready for use in less than 1 hour.
- Easier to use than other 200 or 300 hole systems.
- Dishwasher safe Inox (SS316) powder contact parts. Filler base in lightweight, corrosion resistant anodized aluminum.
- Change capsule sizes in only 5 minutes. Change Parts for Size 000-4 & 00el-2el.
- Orienter Base accommodates capsule sizes 00-4. Size 4 requires Bottom Plate 7180-01 for usage.
- All metal parts Inox (SS316) for reduced static
- Adjustable capsule body height to make over-filling easier & faster. Unique features include lightweight 2-piece orienter (3700), orienter locator and oversize Powder Tray & Push Plate.
- Designed for fast disassembly & visual verification of cleaning.
- Adjustable capsule body height to make over-filling easier & faster.
- Fill weights, USP/EP method within  $\pm 1-3\%$  of average. Typically all 300 capsules/batch are within  $\pm 4-7\%$  of target. More info at <http://www.profiller.com/3fw>.
- Works with all capsules including gelatin, vegetarian/HPMC and starch capsules.
- Tamper size range optimized for higher fill weights and faster filling.
- Indicator gauges to assist in capsule orientation, separation and locking.
- 3700 fully locks capsules after filling, unlike most hand-held fillers - saving you time and frustration!
- Video tips at [www.profiller.com](http://www.profiller.com). Well-illustrated manuals in multiple languages.
- Compliance: CE, cGMP (validation package available upon request)

**ADDITIONAL INFORMATION:**

**Warranty:** 1 year

**Applicable Products:**

- 📄 7540: Partial Batch Kit, ProFiller 3600/3700 (Sizes 00-4)
- 📄 7599: Vibrating Table, ProFiller 3700 (110-120V)

**Description of Features:**

**ProFiller 3700 & 3800 Filler**

- Ease of Use
- ◆ GMP/Cleaning
- ▲ Speed

◆	SS316 (Inox) Locking Plate. No rubber pad to wear out or clean
●	Rear locking handle for fast, easy movement across front of machine. New internal gear design for easier locking
◆	Dishwasher safe Inox (SS316) powder contact parts including Bottom Plate (not visible)
▲	Adjustable capsule body height for easier & faster over-filling of capsules
◆	Base Frame in light-weight anodized aluminum or upgrade to Inox (SS316)
●	Cam position indicator, great for training new operators
●	Over-size thumb screws for easy size change over
●	Dual springs make Tamper easier to lower & raise
▲	Tamper pin size range optimized for faster filling & higher fill weights. Size ranges: 000 / 00 / 0-1 / 2-3 / 4 / 5

- ◆ Body Sheets fully separate for easy cleaning & visual verification of cleaning
- ◆ Cover for gear mechanism & Vibrator to minimize powder entry (inside view)
- ▲ Powder Tray is double-depth for faster Spreading
- ▲ Multi-use Push Plate speeds up operations

**3700 Orienter**

**3800 Orienter**

●	2-piece Hand-Held Orienter. Shake only the lightweight Orienter Tray	●	Orienter Base common for Size 00-4. Traditional European style Orienter for operators who prefer it
●	Locking indicator confirms capsules are fully locked	●	Guide lines to locate Caps Tray for partial batches
▲	Locating holes in Orienter and tapered locating pins for fast & accurate orientation of capsules into Caps Tray	●	Locking indicator confirms capsules are fully locked

**Adjustment-Free  
Size Changeovers**

**Shipped Fully Assembled**

**Dimensions (For Powder Hood Sizing):**

