



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8110




DRY BLENDING CONTAINER w/LID (500 mL)
Model SS115, Stainless Steel

DESCRIPTION: This container is required with the Waring Blender Base, to achieve the various dry blending requirements.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	500
Material (Lid & Container)	Stainless Steel

REQUIRED PRODUCT:

-  8106: Blender Base (1L / 1.2 L)