



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8116



BEAKER, GRADUATED GLASS (4000 mL, Heavy Duty Glass)

Graduation Range: 500–3500 mL, Graduation Interval: 250 mL, Dual Scale

DESCRIPTION: This heavy duty beaker with spout is specifically designed to offer good mechanical strength under harsh conditions, such as mechanized washing operations. For convenience, it has a double graduated metric scale to indicate approximate content and has an extra-large marking spot. It also conforms to ASTM E-960 standards.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	4000
Material	Type I Borosilicate glass
Graduation Range (mL)	500 – 3500
Graduation Interval (mL)	250
Accuracy	± 5% of full capacity
O.D. × H (mm)	160 × 250