



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING





PRODUCT NO.: 8260



Balance Accessories, Calibration Weight, 2Kg

DESCRIPTION: Primarily used to test commercial weighing devices by state and local weights and measures officials, device installers and service technicians. Class F weights may be used to test most accuracy Class III scales, all scales of Class IIIIL or IIII, and scales not marked with a class designation.

FEATURES:

-  Tolerance: 0.20g
-  Material: 303 Stainless Steel
-  Density: 7.85 g/cm³ at 20°C
-  Dimensions: 2.475" x 3.336"

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM