



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 8263-01**








## **CALIBRATION WEIGHT**

5kg, Cast Iron, Class F w/No Certificate

**DESCRIPTION:** Primarily used to test commercial weighing devices by state and local weights and measures officials, device installers and service technicians. Class F weights may be used to test most accuracy Class III scales, all scales of Class IIIIL or IIII, and scales not marked with a class designation.

### **TECHNICAL SPECIFICATIONS:**

-  Tolerance: 0.50g
-  Style: Grip handle weight
-  Material: 303 Stainless Steel
-  Density: 7.2 g/cm<sup>3</sup> at 20°C
-  Dimensions: 3.25" x 5.38" x 3.88"

### **ADDITIONAL INFORMATION:**