



PRODUCT NO.: 8289



SNAP-ON STOPPER (13mm, Red) Natural Rubber

DESCRIPTION: This design provides maximum sealing for vials with a blow-back glass finish, which is a small ridge on the interior of the vial's neck. The snap-on stopper locks into the blow-back ridge, providing a tight seal and keeping the stopper from backing out. This style is not recommended for vials with a straight-walled finish. The blow-back feature is typically found in vials from North American manufacturers.

FEATURES:

- Contains latex due to natural rubber formulation
- Manufactured from the highest quality raw materials
- Components manufactured to yield exceptional dimensional stability
- Snap-on style provides maximum sealing for vials with a blow-back glass finish
- Autoclavable

TECHNICAL SPECIFICATIONS:

- Durometer: 40
- Style: Snap-On
- Rubber Material: Natural Rubber