



PRODUCT NO.: 8362



DIGITAL POCKET THERMOMETER (-40 TO 150°C / -40 TO 302°F)

5" x 0.138" Ø, Dual Scale, Division: 0.1°C / 0.1°F, 1.5V Oxide Battery Included

DESCRIPTION: These stainless steel pocket thermometers are ideal for use in cuvettes, incubators, soil, petroleum products, foods, laboratory and wherever a quick and accurate measurement of temperature is required.

FEATURES:

- Dimensions, W x D x H (in / cm): 35 x 27 x 1 / 89 x 68.6 x 2.5
- Dual Scale: -40 to 150°C / -40 to 302°F
- °F/°C Switch
- Fast 1-second update time
- Swiveling head that can be moved 180° for easy reading in any position
- ON/OFF switch and the Scale Selection switch are located on the face for simple access
- Comes complete with protective pocket case with clip and 1.5V silver oxide battery
- Automatic power-off after 10 minutes of idle use
- Division: 0.1°C / 0.1°F
- Stem is 5" long x 0.138" diameter
- Accuracy is ±2°F (±1°C)