



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8397



Calibration Weight **400g, Stainless steel, Polycarbonate Case & Glove**

DESCRIPTION: Can be used as a reference standard in calibrating other weights and is appropriate for calibrating high-precision analytical balances with a readability as low as 0.1 mg to 0.01 mg.

FEATURES:

- ✓ Tolerance: 0.975mg
- ✓ Material: 316 Stainless Steel
- ✓ Density: 7.95 g/cm³ at 20°C
- ✓ Dimensions (Weight Only): 1.449" x 1.896"
- ✓ Dimensions (Weight in Case): 3.5" x 4.5"
- ✓ Case: Supplied in high quality polycarbonate case

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM