



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 8437-01**



## **SPARE PARTS KIT (Impak Corp)** For STS-7 Tube Sealer, Single Element


**DESCRIPTION:** This single element spare parts kit is used for replacement of heating elements and Teflon covers after prolonged use.

### **ADDITIONAL INFORMATION:**

#### **Included Parts:**

-  1 x Heating Element 1/4" Wide x 7" Long
-  1 x Precut 7" Seal Rubber
-  1 x Cold Side Teflon Cover
-  1 x Hot Side Teflon Cover
-  1 x 10 mm Thick Teflon Cover

#### **Applicable Products:**

-  9567-01: STS-7 Tube Sealer, Single Element