



PRODUCT NO.: 8684



BEAKER, GRADUATED GLASS TYPE I (10 mL)
No Graduation Range/Interval, Single Scale

DESCRIPTION: This beaker meets ASTM E-960 Type I specifications for classification, design, capacity and dimensions and markings.

FEATURES:

- ✓ Easy-pour-sprout
- ✓ White, permanent graduations
- ✓ Large marking square
- ✓ Made from high quality borosilicate glass with a low coefficient of expansion
- ✓ Manufactured with uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength

TECHNICAL SPECIFICATIONS:

Capacity (mL)	10
Dimensions, OD x Height (mm)	25 x 35
Accuracy	± 5% capacity