



PRODUCT NO.: 8688



SERUM VIAL w/STOPPER & SEAL, STERILE (Amber, 2mL, 13mm)

Glass Type 1, w/Silver Seal

DESCRIPTION: All vials are sterile and pyrogen-free. Prior to vial assembly, vials are dry-heat sterilized and stoppers are steam-sterilized. The vials/stoppers are assembled on an automated line located in a Class 100 suite. The vials are subsequently capped in a Class 100,000 environment. Lot release sterility testing conforms to requirements set forth in USP <71>.

FEATURES:

- ☑ Sterile
- ☑ Volume: 2 mL
- ☑ Neck thread: 13 mm
- ☑ Cap Color: Silver
- ☑ Glass: Amber Type 1
- ☑ Stopper: 4401/45, Synthetic Chlorobutyl, Latex-Free
- ☑ Seal: Open-Top