



HOW IT WORKS

1. Uses tuned leverage to help push any size EMP jar plungers.
2. Simple operation provides efficient use of existing pharmacy materials to keep costs of system and materials lower.

EQUIPMENT & DEVICES

TOPI-CLICK® FILLING STATION™



FEATURES & BENEFITS

- Most EMP jar sizes can be quickly interchanged
- Optimizes the closed EMP system; less batch transfers and minimal storage costs
- Works with right or left handed operators via switchable operating lever
- Balanced actuation arm for easy operation
- Made of durable stainless and powder coated steel for years of operation
- Easy to clean and wipe down
- Contrast background to help assure quality filling in see-through containers
- Stabilizer plate for use when filling adapter type dispensers with EMP jars
- Includes NEW nozzle specific for filling Topi-CLICK®

MEDISCA DISCLAIMER: All information presented here, including specifications and suggested uses, are those of the manufacturer of each specific item/device/product. Any representations, warranties or liability for the products are solely those of the manufacturer. MEDISCA makes no representations and takes no responsibility for the functioning or suitability, or clinical or therapeutic use of the item/device/product with or without active ingredients. The pharmacist is hereby advised to use this product in accordance with a licensed prescriber's prescription. MEDISCA cannot be held liable to any person or entity concerning claims, loss, or damage caused by, or alleged to be caused by the use or misuse of the information contained in this suggested product information sheet.

TOPI-CLICK® DISCLAIMER: The Topi-CLICK® product may be improved over time based on customer feedback on usage. Patent Pending. The Topi-CLICK® Filling Station™ Custom Rx TDA, LLC and/or the patent holders are not liable for any damages arising from the products dispensed in the Topi-CLICK® Filling Station™ or its improper use. References: Topi-CLICK® is a registered trademark of Custom Rx TDA, LLC, Atlanta. **The above notes are for certain batches of the referred to bases and without other ingredients. All bases/batches/ formulas have slight differences and therefore need independent calibration. Some viscous solutions are not practical for use in the Topi-CLICK®.



Also available in classic pink, blue, white and black

EQUIPMENT & DEVICES

TOPI-CLICK[®]



FEATURES & BENEFITS

- Domed applicator pad eliminates the need to use hands to apply topicals
- Audible, visual and tactile click makes it easy for all patients to accurately dose medications
- "Refill" indication for prompting on-time re-buys that help eliminate missed applications
- Lower cost with simplified labeling, which is fast and easy as a prescription bottle
- Easier to fill and use with more control over plunger movement than a topical syringe
- Eight colors for differentiating medications (also available in pink, blue, white and new black)
- Equipped with a flip tab for easier leverage to pop the top
- Applicable for total dispense weight of approximately 5g to 35g of topical preparations
- Turn the base "1 click" (a ¼ turn) to dispense one metered amount of 0.25mL
- Made of FDA medical-grade polypropylene material
- UV-block plastic protects contents from UV rays

CALL NOW OR VISIT OUR WEBSITE TO PLACE AN ORDER

 **1.800.932.1039** | **WWW.MEDISCA.COM** | **1.800.665.6334** 