



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 8719**










**SERUM VIAL w/STOPPER & SEAL, STERILE (Amber, 30mL, 20mm)**

Glass Type 1, w/Silver Seal

**DESCRIPTION:** All vials are sterile and pyrogen-free. Prior to vial assembly, vials are dry-heat sterilized and stoppers are steam-sterilized. The vials/stoppers are assembled on an automated line located in a Class 100 suite. The vials are subsequently capped in a Class 100,000 environment. Lot release sterility testing conforms to requirements set forth in USP <71>.

**FEATURES:**

-  Sterile
-  Volume: 30 mL
-  Neck thread: 20 mm
-  Cap Color: Silver
-  Glass: Amber Type 1
-  Stopper: 4401/45, Synthetic Chlorobutyl, Latex-Free
-  Seal: Open-Top