



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8720




Vial, Amber Glass 10mL, 20mm (Sterile)

DESCRIPTION: Prior to vial assembly, vials are dry-heat sterilized and stoppers are steam-sterilized. The vials/stoppers are assembled on an automated line located in a Class 100 suite. The vials are subsequently capped in a Class 100,000 environment. Lot release sterility testing conforms to requirements set forth in USP <71>.

TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	12.5 ± 0.5
Height (mm)	52.0 ± 0.5
Thickness of bottle-lip (mm)	3.9 ± 0.2
Height of bottle-neck (mm)	3.5 ± 0.3
Effective height of bottle-body (mm)	40.0 ± 1.0
Outer diameter of bottle-mouth (mm)	19.8 ± 0.2
Inner diameter of bottle-mouth (mm)	12.6 ± 0.2
Outer diameter of bottle-neck (mm)	15.8 ± 0.3
Body diameter (mm)	22.0 ± 0.2
Material	Borosilicate glass



USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM